































# PROGRAMMING CABLES AND ADAPTERS





















	<b>RITE</b>	<b>Micro RITE</b>	<b>Nano RITE</b>	<b>Pico RITE</b>				
				 390-01-180-05				
	<b>BTEs</b>	<b>Nano</b>			<b>Micro / Compact Power/Super Power</b>			<b>Power BTE</b>
			 390-01-151-00			 399-50-640-00		
<p><b>LEFT:</b> HiPro 384-20-033-00* NOAH Link: 384-20-035-00</p> <p><b>RIGHT:</b> HiPro 384-20-032-00* NOAH Link: 384-20-034-00</p>	<b>ITCD Chronos Inizia</b>	 Flex Connect with battery pill 390-01-040-00						
 <p>FittingLINK All wireless BTE and Custom Hearing Instruments</p>		<b>Chronos Inizia</b> All Custom Hearing Instruments	 Flex Connect 390-01-180-05					
	<b>Saphira Juna Carista Acuva</b> All Custom Hearing Instruments	 Mini Flex Connect 117468						

**Manufacturer:**  
Bernafon AG  
Morgenstrasse 131  
3018 Bern  
Switzerland  
www.bernafon.com

**Local Manufacturer  
& Distributor:**  
Bernafon Canada  
500 Trillium Drive, Unit 15  
Kitchener, ON, N2R 1A7  
www.bernafon.ca

# ACOUSTIC OPTIONS OVERVIEW

JUNA 9 | 7  
ACRIVA 9 | 7  
SAPHIRA 5 | 3  
CARISTA 5 | 3  
INIZIA 3 | 1  
NEVARA 1

				Instant Fitting							Custom Fitting			
				Open Dome	Bass Dome, Double Vent	Bass Dome, Single Vent	Power Dome	Open Dome	Tulip	Power Dome	Micro/Canal Mold	Lite Tip	Power Mold	
														
BTE	Standard Earhook				N/A	N/A	N/A	N/A	N/A	N/A	Standard earmolds as recommended by clinician.			
	miniFit Thin Tube 0.9/1.3 mm*				6 mm 8 mm 10 mm	6 mm 8 mm 10 mm 12 mm	6 mm 8 mm 10 mm 12 mm	6 mm 8 mm 10 mm 12 mm	N/A	N/A	N/A	Occluded 1.0 mm 1.4 mm 2.4 mm 3.0 mm	Occluded 0.7 mm (long) 0.7 mm (short) 1.2 mm 1.5 mm	N/A
RITE	60, 85 Speaker				6 mm 8 mm 10 mm	6 mm 8 mm 10 mm 12 mm	6 mm 8 mm 10 mm 12 mm	6 mm 8 mm 10 mm 12 mm	N/A	N/A	N/A	Occluded 1.0 mm 1.4 mm 2.4 mm 3.0 mm	Occluded 0.7 mm (long) 0.7 mm (short) 1.2 mm 1.5 mm	N/A
	100 Speaker				N/A	6 mm 8 mm 10 mm 12 mm	6 mm 8 mm 10 mm 12 mm	6 mm 8 mm 10 mm 12 mm	N/A	N/A	N/A	Available on request dependent on ear size	N/A	Occluded 0.8 mm 1.4 mm 2.4 mm 3.0 mm Other **
	M-Speaker bent				N/A	N/A	N/A	N/A	6 mm 8 mm 10 mm	12 mm	6 mm 8 mm 10 mm 12 mm	Occluded 1.0 mm 1.4 mm 2.4 mm 3.0 mm	Occluded 0.7 mm (long) 0.7 mm (short) 1.2 mm 1.5 mm	N/A
	P-Speaker				N/A	N/A	N/A	N/A	N/A	12 mm	6 mm 8 mm 10 mm 12 mm	N/A	N/A	Occluded 1.0 mm 1.4 mm 2.4 mm

\* Spira Flex thin tubes must be used with a BTE adapter which replaces the standard earhook.

\*\* Please note: Instant fit domes are specific for spira flex vs. M-Speaker vs. Pico mini speakers. These domes ARE NOT interchangeable.